

SOP-103: DART Editing for Bacterial Genome Engineering via Conjugation

Organism: ____

Version: 1.0

Effective Date: ____

Author:

Approved by: ____

1. PURPOSE AND SCOPE

1.1 Purpose

Project-specific purpose:

1.2 Process Overview

This SOP describes RNA-guided DNA insertion using CRISPR-Associated Transposase (CAST) systems delivered via bacterial conjugation for genome engineering of non-model gut microbiome bacteria.

CAST-Mediated Integration Workflow:

PREREQUISITE: Vector Construction - gRNA design and cloning completed via Gibson Assembly SOP - DART vector verified and transformed into donor strain - Donor culture prepared and ready for conjugation - **Note:** This SOP begins with a verified DART plasmid in donor cells

STEP 1: Conjugation-Based Delivery - Prepare recipient bacterial culture from gut microbiome - Perform bacterial mating between donor and recipient cells - Plate on selective medium (no DAP, with integration marker) - Incubate at 30°C for 24-48 hours - **Purpose:** Transfer CAST machinery to target bacteria

STEP 2: Selection and Clonal Isolation - Select transconjugants on antibiotic plates - Screen colonies by external-external PCR - Re-streak positive colonies to ensure clonality - Isolate clonal integrants with only integration allele - **Purpose:** Obtain pure clones with genomic insertions

STEP 3: Integration Verification - PCR amplification of genome-transposon junctions - Sanger sequencing to confirm integration site and orientation - Verify 5-bp target site duplication (TSD) - Confirm absence of wild-type allele - **Purpose:** Validate on-target integration events

Key Advantages of CAST Systems: - RNA-programmable targeting (no fixed recognition sequences) - High integration efficiency (50-99% for 1-2 kb payloads) - No requirement for homologous recombination machinery - Suicide vector eliminates need for plasmid curing - Compatible with diverse Gram-negative bacteria - Stable, single-copy chromosomal integration

Mechanism: CAST systems use a TniQ-Cascade complex to bind 32-bp target sequences via crRNA guide. The TnsABC transposase catalyzes non-replicative cut-and-paste insertion of the mini-transposon ~49 bp downstream of the target site. Integration generates 5-bp target site duplications (TSDs) and predominantly T-RL orientation products (>90%).

2. PROJECT INFORMATION

Project Name: ____

Target Bacterium/Bacteria: - Source: ____ - *Species/Strain (if known):* ____ - Gram stain: ☐ Negative ☐ Positive ☐ Unknown

Genomic Integration Site: ____ - *Target gene/locus:* ____ - *Genomic coordinates:* ____ **to** ____ - Integration purpose: ☐ Gene knockout ☐ Reporter insertion ☐ Promoter capture

DART Vector System: ____ - *Backbone:* ____ - *CAST system:* ____ - *Plasmid selection:* ____ - *Integration selection:* ____

Payload Information: - Size: ____ **bp** - **Components:** ____ - Expected orientation: ☐ T-RL (preferred) ☐ T-LR ☐ Either

Donor Strain: ____ - *Strain:* ____ - *Genotype:* ____ - DAP requirement: ☐ Yes ☐ No

Recipient Strain(s): - Strain ID: ____ - *Isolation source:* ____ - *Growth temperature:* 30°C - *Growth medium:* ____

3. MATERIALS AND EQUIPMENT

3.1 Reagents and Media

Category	Item	Specification / Notes
Plasmids & Strains	<input type="checkbox"/> DART plasmid construct	Vector with custom gRNA and payload
	<input type="checkbox"/> Donor strain competent cells	DAP-dependent (e.g., WM3064, MFDpir+)
	<input type="checkbox"/> Recipient bacterial culture	Target bacteria for genome editing
Antibiotics (1000× stocks)	<input type="checkbox"/> Spectinomycin (100 mg/mL)	For donor selection during cloning
	<input type="checkbox"/> Selection antibiotic	As appropriate for transconjugant selection
	<input type="checkbox"/> Diaminopimelic acid (DAP, 50 mg/mL)	For donor strain growth (required)
Growth Media	<input type="checkbox"/> LB broth	Standard lysogeny broth
	<input type="checkbox"/> LB agar plates	For conjugation and selection
	<input type="checkbox"/> LB + DAP (0.3 mM final)	For donor strain maintenance
	<input type="checkbox"/> LB + appropriate antibiotics	For selection as specified
	<input type="checkbox"/> SOC or LB (recovery medium)	For post-transformation recovery
	<input type="checkbox"/> Recipient-specific growth medium	As required for target organism
DNA Purification Kits	<input type="checkbox"/> Plasmid miniprep kit	For plasmid extraction
	<input type="checkbox"/> PCR purification kit	For PCR product cleanup
	<input type="checkbox"/> Genomic DNA extraction kit	For gDNA extraction
General Lab Supplies	<input type="checkbox"/> PBS (1×, pH 7.4)	Gibco #10010023 or make fresh
	<input type="checkbox"/> Sterile 0.22 µm filters	For media sterilization
	<input type="checkbox"/> Sterile spreader or plating beads	For plating bacteria
	<input type="checkbox"/> Pipette tips (sterile)	P2, P10, P20, P200, P1000
	<input type="checkbox"/> Microcentrifuge tubes (1.7 mL)	For reactions and storage
	<input type="checkbox"/> 50 mL conical tubes	For cultures
	<input type="checkbox"/> Petri dishes (100 × 15 mm)	Sterile, for plating
	<input type="checkbox"/> Cryovials	For glycerol stocks
	<input type="checkbox"/> Sterile inoculation loops	Single-use plastic loops

4. SAFETY AND PRECAUTIONS

- ☐ Wear appropriate PPE (lab coat, gloves, safety glasses)
- ☐ Work in biosafety cabinet when handling bacterial cultures
- ☐ Mealworm gut microbiome may contain BSL-2 organisms—follow institutional biosafety protocols
- ☐ UV protection when using gel transilluminator
- ☐ Gentamycin is toxic—handle with care, avoid inhalation or skin contact

- [] Autoclave all biological waste before disposal
 - [] Follow institutional guidelines for GMO handling and disposal
 - [] DAP-dependent donor strain cannot survive outside laboratory—contains biocontainment feature
-

5. PART 1: CONJUGATION-BASED DELIVERY

7.1 Overview

Transfer pBFC0619-gRNA from *E. coli* WM3064 donor to recipient gut bacteria via bacterial conjugation. Select for transconjugants on gentamycin (no DAP), preventing donor growth.

Mating Method: ☐ Filter mating ☐ Plate mating (direct mix)

7.2 Prepare Recipient Bacterial Culture

Date: ____

Recipient Strain: ____

1. Inoculate fresh recipient culture from frozen stock or plate
2. Grow in ____ mL medium
3. Temperature: 30°C
4. Time: ____ hours (to mid-log phase, OD600 = 0.4-0.8)
5. Measure OD600: ____

Growth Conditions: - Medium: ____ - Incubation time: ____ hours - Starting OD600: ____ - Final OD600: ____ (mid-log preferred)

✓ **QC Checkpoint:** Recipient culture in exponential growth phase (OD600 = 0.4-0.8)

5.3 Perform Bacterial Conjugation (Plate Mating Method)

Date: ____

Mating Reaction:

1. Prepare donor cells:

- Harvest 1 mL donor culture
- Pellet at 5,000g for 3 minutes
- Wash with 1 mL sterile PBS (pH 7.4)
- Pellet again, resuspend in 100 µL PBS
- Measure OD600: ____

2. Prepare recipient cells:

- Harvest 1 mL recipient culture (from Step 5.2)
- Pellet at 5,000g for 3 minutes
- Wash with 1 mL sterile PBS
- Pellet again, resuspend in 100 µL PBS
- Measure OD600: ____

3. Prepare mating mixtures:

Mix donor and recipient (Conjugation):

- Combine donor and recipient in 1:1 ratio (by OD or volume)
- Example: 50 µL donor + 50 µL recipient
- Mix gently by pipetting

Donor only (Control):

- 50 µL donor cells in PBS

Recipient only (Control):

- 50 µL recipient cells in PBS

4. Spot on conjugation plate:

- Prepare LB agar plate (no antibiotics, with DAP if donor requires)
- Spot each mixture separately on the same plate (label positions):
 - **Conjugation mix:** 10-50 µL in position 1
 - **Donor only:** 10-50 µL in position 2
 - **Recipient only:** 10-50 µL in position 3
- Allow all spots to dry completely in biosafety cabinet (10-15 min)

5. Incubate for mating:

- Incubate plate at 30°C for 6-24 hours
- Optimal time: ___ hours (typically 18-24 h)

Conjugation Conditions: - Donor:Recipient ratio: ___ : - Mating time: ___ hours - Temperature: 30°C

✓ **QC Checkpoint:** All three mating spots visible on plate after incubation

5.4 Select for Transconjugants

Date: ___

1. Harvest mating spots:

- **Conjugation spot:** Add 500 µL sterile PBS, resuspend by pipetting, transfer to tube labeled "Conjugation"
- **Donor-only spot:** Add 500 µL sterile PBS, resuspend by pipetting, transfer to tube labeled "Donor only"
- **Recipient-only spot:** Add 500 µL sterile PBS, resuspend by pipetting, transfer to tube labeled "Recipient only"

2. Prepare serial dilutions for each sample:

- 10^0 : Use 100 µL undiluted
- 10^{-1} : 100 µL sample + 900 µL PBS
- 10^{-2} : 100 µL from 10^{-1} + 900 µL PBS
- 10^{-3} : 100 µL from 10^{-2} + 900 µL PBS

3. Plate on selective medium:

For **each sample** (Conjugation, Donor-only, Recipient-only), plate 100 µL of each dilution on:

- **Selection plates:** LB + ___ (antibiotic), NO DAP ← transconjugants only
- **Total count plates (optional):** LB only (no antibiotic) ← total viable cells

Spread evenly with sterile beads or spreader. Allow to dry completely.

4. Expected results:

- **Conjugation plates (LB + antibiotic, no DAP):** Transconjugant colonies
- **Donor-only plates (LB + antibiotic, no DAP):** NO growth (donor requires DAP)
- **Recipient-only plates (LB + antibiotic, no DAP):** NO growth (no resistance marker)

5. Incubate:

- Temperature: 30°C
- Time: 24-48 hours (until colonies visible)

Selection Conditions: - Antibiotic: ___ at µg/mL - DAP included: ☐ Yes ☐ No - Incubation temperature: 30°C - Incubation time: ___ hours

Colony Counts:

Sample	Plate Type	Dilution	Colonies Counted	CFU/mL	Notes
Conjugation	LB + antibiotic (no DAP)	10 ⁰	___	___	Transconjugants
Conjugation	LB + antibiotic (no DAP)	10 ⁻¹	___	___	Transconjugants
Conjugation	LB + antibiotic (no DAP)	10 ⁻²	___	___	Transconjugants
Conjugation	LB only	10 ⁰	___	___	Total cells
Donor only	LB + antibiotic (no DAP)	All	___	___	Should be 0
Recipient only	LB + antibiotic (no DAP)	All	___	___	Should be 0

Conjugation Efficiency: ___ % (transconjugants / total cells × 100)

✓ **QC Checkpoint:** Transconjugant colonies visible on conjugation plates, NO growth on donor-only and recipient-only control plates

Note: If donor-only plates show growth, DAP may have contaminated the medium. If recipient-only plates show growth, the recipient may have native antibiotic resistance—test a higher antibiotic concentration.

6. PART 2: COLONY SELECTION AND CLONAL ISOLATION

6.1 Overview

Screen transconjugants by colony PCR to identify clones with successful genomic integration. Re-streak to ensure clonality, as initial colonies may be heterogeneous.

Expected Integration Products: - T-RL orientation (preferred, >90%): Transposon right end proximal to target - T-LR orientation (rare, <10%): Transposon left end proximal to target - 5-bp target site duplication (TSD) flanking insertion

6.2 Screen Colonies by PCR

Date: ___

Screening Strategy: External-external PCR (genome primers flanking insertion site)

1. Pick colonies:

- Select 10-20 well-isolated colonies from LB + antibiotic plates
- Resuspend each in 40 µL sterile water in PCR strip tubes
- Label tubes: Colony #1, #2, #3, etc.

2. Archive colonies:

- Using fresh LB + antibiotic plate, spot 1 µL from each resuspension
- Allow to dry, then incubate at 30°C overnight as backup

3. Prepare cell lysate:

- Heat resuspensions at 95°C for 10 minutes
- Cool to room temperature
- Dilute 1:20 in water (2 µL lysate + 38 µL water)

Colony Lysates:

Colony #	Lysate Prepared?	Spotted on Backup Plate?
1-5	<input type="checkbox"/>	<input type="checkbox"/>
6-10	<input type="checkbox"/>	<input type="checkbox"/>
11-15	<input type="checkbox"/>	<input type="checkbox"/>
16-20	<input type="checkbox"/>	<input type="checkbox"/>

6.3 External-External PCR for Integration Detection

Date: __

Primer Design: - **UpstreamFwd:** Binds genome ~500 bp upstream of target site - **DownstreamRev:** Binds genome ~500 bp downstream of integration site - Expected WT product: ~__ bp (no insertion) - **Expected integration product:** ~ bp (includes payload + flanking sequence)

Primers:

Primer Name	Sequence (5' → 3')	Tm	Location
Upstream_Fwd	__	__ °C	Genome, upstream
Downstream_Rev	__	__ °C	Genome, downstream

PCR Reaction (25 µL per colony):

Component	Volume per Reaction	Mastermix for 22 Reactions
Diluted lysate	2 µL	N/A
Upstream_Fwd (10 µM)	1.25 µL	27.5 µL
Downstream_Rev (10 µM)	1.25 µL	27.5 µL
Q5 Reaction Buffer (5×)	5 µL	110 µL
dNTPs (10 mM)	0.5 µL	11 µL
Q5 DNA Polymerase	0.25 µL	5.5 µL
Nuclease-free water	14.75 µL	324.5 µL

Prepare mastermix for all reactions plus 2 extra (20 colonies + WT control + no-template control)

PCR Cycling Conditions:

1. 98°C - 30 sec (initial denaturation)
2. 30 cycles:
 - 98°C - 10 sec
 - ____ °C - 20 sec (annealing, adjust based on primer Tm)
 - 72°C - 120 sec (extension)
3. 72°C - 2 min (final extension)
4. 4°C - hold

Controls: - **WT control:** Genomic DNA from recipient strain (no integration) - Expected band: ~__ bp (no insertion) - **No-template control:** Water instead of lysate - Expected: No band

6.4 Gel Electrophoresis and Analysis

Date: __

1. Add 5 µL of 6× loading dye to each PCR reaction

2. Load 10 µL per well on 1% agarose gel (0.8% for larger products)
3. Include 1 kb DNA ladder in first and last lanes
4. Run at 120 V for 40-50 minutes (until bands well-separated)
5. Image gel using UV or blue-light transilluminator

Expected Results: - **Integration positive:** Single band at ~__ bp (includes payload) - **Wild-type:** Band at ~__ bp (no insertion) - **Heterogeneous:** Both WT and integration bands present - **No integration:** Only WT band

Screening Results:

Colony #	WT Band?	Integration Band?	Interpretation	Select for Re-streak?
1	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Integration <input type="checkbox"/> WT <input type="checkbox"/> Mixed	<input type="checkbox"/> Yes <input type="checkbox"/> No
2	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Integration <input type="checkbox"/> WT <input type="checkbox"/> Mixed	<input type="checkbox"/> Yes <input type="checkbox"/> No
3	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Integration <input type="checkbox"/> WT <input type="checkbox"/> Mixed	<input type="checkbox"/> Yes <input type="checkbox"/> No
4	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Integration <input type="checkbox"/> WT <input type="checkbox"/> Mixed	<input type="checkbox"/> Yes <input type="checkbox"/> No
5	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Integration <input type="checkbox"/> WT <input type="checkbox"/> Mixed	<input type="checkbox"/> Yes <input type="checkbox"/> No
6-20

✓ **QC Checkpoint:** At least 3-5 colonies show integration band

6.5 Re-streak for Clonal Purity

Date: __

Rationale: Initial colonies are often heterogeneous (contain both integrated and WT alleles). Re-streaking ensures clonal populations.

1. Select 3-5 colonies showing integration bands (from Step 6.4)
2. From backup plates, pick well-isolated colonies
3. Streak for single colonies on fresh LB + antibiotic plates
4. Incubate at 30°C for 24-48 hours

Re-streaked Clones:

Original Colony #	Re-streak Plate ID	Incubation Time
Colony #__	Plate A	__ hours
Colony #__	Plate B	__ hours
Colony #__	Plate C	__ hours

1. After growth, pick 3-5 single colonies from each re-streak plate
2. Repeat colony PCR screening (Steps 6.2-6.4)
3. Identify clones with ONLY integration band (no WT band)

Second-Round Screening Results:

Re-streak Plate	Colony #	WT Band?	Integration Band?	Clonal?
Plate A	1	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No
Plate A	2	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No
Plate B	1	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No
Plate B	2	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No
Plate C	1	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No

✓ **QC Checkpoint:** At least one clonal integrant identified (integration band only, no WT)

7. PART 3: INTEGRATION VERIFICATION

7.1 Overview

Confirm on-target integration by sequencing genome-transposon junctions. Verify integration orientation (T-RL vs T-LR) and precise insertion site.

7.2 Junction PCR for Sequencing

Date: ____

Strategy: Amplify genome-transposon junctions using genome-specific and transposon-specific primers

Primer Pairs:

5' Junction (Upstream): - *GenomeFwd*: Binds genome ~300 bp upstream of target - *TnsRRev*: Binds transposon right end

3' Junction (Downstream): - *TnsLFwd*: Binds transposon left end
- *GenomeRev*: Binds genome ~300 bp downstream

Primers:

Primer Name	Sequence (5' → 3')	Tm	Target
Genome_Fwd	____	__ °C	Genome
TnsR_Rev	5'-CGCCAGGGTTTTCCCAGTC-3'	60°C	Transposon R end
TnsL_Fwd	5'-CTGTAGCGCGTTTTTCATCGG-3'	60°C	Transposon L end
Genome_Rev	____	__ °C	Genome

Select 2-3 clonal integrants from Step 6.5

1. Inoculate clones in 5 mL ____ medium + antibiotic
2. Grow overnight at 30°C
3. Extract genomic DNA using Wizard Genomic DNA Purification Kit
4. Quantify DNA

gDNA Yields:

Clone ID	Concentration (ng/μL)	A260/A280	Pass/Fail
Clone # ____	____	____	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
Clone # ____	____	____	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
Clone # ____	____	____	<input type="checkbox"/> Pass <input type="checkbox"/> Fail

PCR Reaction for Each Junction (25 µL):

Component	Volume
gDNA (10-50 ng)	2 µL
Forward primer (10 µM)	1.25 µL
Reverse primer (10 µM)	1.25 µL
Q5 Reaction Buffer (5×)	5 µL
dNTPs (10 mM)	0.5 µL
Q5 DNA Polymerase	0.25 µL
Nuclease-free water	14.75 µL

PCR Cycling Conditions:

1. 98°C - 30 sec
2. 30 cycles:
 - 98°C - 10 sec
 - ____°C - 20 sec (annealing)
 - 72°C - 30 sec (extension: ~300-500 bp product)
3. 72°C - 2 min
4. 4°C - hold

Run on 1.5% agarose gel - Expected 5' junction product: ~__ bp - Expected 3' junction product: ~__ bp

7.3 Gel Extract and Sequence Junctions

Date: __

1. Excise junction PCR products from gel
2. Gel extract using QIAquick Gel Extraction Kit
3. Elute in 30 µL elution buffer
4. Quantify DNA

Purified Junction Products:

Clone ID	Junction	Concentration (ng/µL)	Pass/Fail
Clone #__	5'	__	<input type="checkbox"/> Pass <input type="checkbox"/> Fail
Clone #__	3'	__	<input type="checkbox"/> Pass <input type="checkbox"/> Fail

1. Submit for Sanger sequencing
 - Use PCR primers (*GenomeFwd*, *GenomeRev*) for sequencing
 - Sequence both junctions for each clone

7.4 Analyze Sequencing Results

Date: __

Expected Features: - **5-bp Target Site Duplication (TSD):** Genomic sequence duplicated on both sides of insertion - **Integration distance:** ~49 bp downstream of target site (3' end of 32-bp target) -

Orientation: T-RL (right end proximal) or T-LR (left end proximal) - **Intact transposon ends:** Proper junction between genome and transposon

Sequencing Analysis:

Clone ID	5' Junction Correct?	3' Junction Correct?	TSD Identified?	Distance from Target	Orientation	Pass/Fail
Clone #__	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	_ bp	<input type="checkbox"/> T-RL <input type="checkbox"/> T-LR	<input type="checkbox"/> PASS <input type="checkbox"/> FAIL
Clone #__	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	_ bp	<input type="checkbox"/> T-RL <input type="checkbox"/> T-LR	<input type="checkbox"/> PASS <input type="checkbox"/> FAIL
Clone #__	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	<input type="checkbox"/> Yes <input type="checkbox"/> No	_ bp	<input type="checkbox"/> T-RL <input type="checkbox"/> T-LR	<input type="checkbox"/> PASS <input type="checkbox"/> FAIL

5-bp TSD Sequence: 5'-__-3'

Validated clone(s): Clone #_, Clone #, Clone #_

Selected for archiving: Clone #__

✓ **Final QC Checkpoint:** At least one fully validated clone with correct junctions and on-target integration

8. STRAIN ARCHIVING

Date: __

Validated Clone ID: Clone #__

Strain Name: __

Genotype: - Parent strain: _ - *Modification:* - *Integration site:* - *Payload:* _ - Orientation: ☐ T-RL ☐ T-LR

Description:

Glycerol Stock Preparation:

1. Inoculate validated clone in 5 mL __ medium + antibiotic
2. Grow overnight at 30°C to saturation (OD600 >2)
3. Mix 500 µL culture + 500 µL sterile 50% glycerol in cryovial
4. Vortex thoroughly to mix
5. Label cryovials clearly
6. Flash freeze in liquid nitrogen or dry ice/ethanol bath
7. Store at -80°C

Storage Locations:

Vial	Box	Position	Freezer	Date
1	__	__	__	__
2	__	__	__	__
3	__	__	__	__

Backup Location: __

Lab Notebook Reference: __

Electronic Files: - Sequencing results: _ - *Gel images:* - *Plasmid maps:* _

9. TROUBLESHOOTING

Problem	Possible Cause	Solution
No colonies after gRNA cloning	<ul style="list-style-type: none"> Incomplete digestion Ligation failed Non-competent cells 	<ul style="list-style-type: none"> Verify complete digestion by gel Check ligase activity with control Use fresh competent cells
High background on no-insert control	<ul style="list-style-type: none"> Vector self-ligation Incomplete digestion 	<ul style="list-style-type: none"> Increase digestion time to 3-4 hours Add alkaline phosphatase (CIP) to digest Use higher insert:vector ratio
Incorrect gRNA sequence	<ul style="list-style-type: none"> Wrong oligos ordered Cloning error 	<ul style="list-style-type: none"> Re-design and re-order oligos Sequence more clones Check for mutations in synthesis
Donor strain grows on -DAP plates	<ul style="list-style-type: none"> DAP contamination Revertant strain 	<ul style="list-style-type: none"> Prepare fresh media and plates Verify donor strain genotype Use fresh donor stock
No transconjugants after conjugation	<ul style="list-style-type: none"> Donor/recipient not viable Insufficient mating time Wrong antibiotic concentration Recipient naturally resistant 	<ul style="list-style-type: none"> Verify donor and recipient OD600 Extend mating time to 24 hours Optimize antibiotic concentration (MIC) Test recipient on antibiotic plates
Donor grows on selection plates (no DAP)	<ul style="list-style-type: none"> DAP contamination in medium Integration into donor genome 	<ul style="list-style-type: none"> Remake selection plates without DAP Verify plates are DAP-free Use fresh donor culture
Low conjugation efficiency (<0.1%)	<ul style="list-style-type: none"> Poor mating conditions Wrong donor:recipient ratio Incompatible strains 	<ul style="list-style-type: none"> Optimize mating time (6-24 h) Try different ratios (1:1, 1:10, 10:1) Increase mating temperature slightly Use filter mating method instead
All colonies are wild-type (no integration)	<ul style="list-style-type: none"> gRNA not functional Transposase not expressed Target site not accessible 	<ul style="list-style-type: none"> Re-design gRNA to different target Verify payload integrity by sequencing Test different genomic sites Extend incubation time at 30°C
Integration efficiency very low (<10%)	<ul style="list-style-type: none"> Suboptimal temperature Large payload Target site chromatin structure 	<ul style="list-style-type: none"> Incubate at 25-30°C instead of 37°C Screen more colonies (30-50) Try alternative target sites Extend incubation to 48 hours
Heterogeneous colonies (mixed WT + integration)	<ul style="list-style-type: none"> Insufficient growth time Clones not pure 	<ul style="list-style-type: none"> Re-streak colonies 2-3 times Pick smaller, younger colonies Extend incubation between re-streaks
PCR gives no product	<ul style="list-style-type: none"> Primers don't bind Poor DNA template quality Wrong annealing temperature 	<ul style="list-style-type: none"> Verify primer sequences Re-extract DNA or use fresh lysate Optimize annealing temperature Use touchdown PCR
Wrong integration site (off-target)		

Problem	Possible Cause	Solution
	<ul style="list-style-type: none"> Off-target sites in genome Non-specific gRNA 	<ul style="list-style-type: none"> Re-design gRNA with BLAST check Avoid targets with >2 off-targets Sequence more clones to find on-target
Both T-RL and T-LR orientations detected	<ul style="list-style-type: none"> Normal behavior for CASTs T-RL is preferred (>90%) 	<ul style="list-style-type: none"> Screen more colonies to find T-RL Both orientations are functional Select preferred orientation if needed
No 5-bp TSD detected in sequence	<ul style="list-style-type: none"> Sequencing quality issue Not a CAST integration event 	<ul style="list-style-type: none"> Re-sequence with higher quality Verify integration by junction PCR Check for contamination
Integration product includes vector backbone	<ul style="list-style-type: none"> Cointegrate formation Incomplete excision 	<ul style="list-style-type: none"> This is rare for Type I-F CASTs Verify by long-read sequencing Use alternative clone
Antibiotic selection doesn't work	<ul style="list-style-type: none"> Recipient has native resistance Antibiotic degraded 	<ul style="list-style-type: none"> Determine MIC for recipient strain Use higher concentration Make fresh antibiotic stocks Try alternative selection marker

10. TIMELINE

Estimated timeline for complete protocol:

Week	Activities	Checkpoints
Week 1	Prepare donor strain with DART plasmid	Verified donor culture ready
Week 2	Grow recipient bacteria, perform conjugation	Transconjugant colonies
Week 3	Screen by PCR, re-streak for clonality	Clonal integrants identified
Week 4	Junction PCR, sequencing, validation	Validated strain archived

Total time: ~4 weeks (from donor preparation to validated strain)

Hands-on time: ~6-8 days total

Note: Timeline assumes: - DART plasmid already constructed and verified (Gibson Assembly SOP completed) - Sanger sequencing turnaround: 1-2 days - Recipient bacteria can be cultured to mid-log in <24 hours

11. REFERENCES

- Gelsinger, D. R. et al. (2024). Bacterial genome engineering using CRISPR-associated transposases. *Nature Protocols*, 19(3), 752-790. <https://doi.org/10.1038/s41596-023-00927-3>
- Klompe, S. E., Vo, P. L. H., Halpin-Healy, T. S. & Sternberg, S. H. (2019). Transposon-encoded CRISPR–Cas systems direct RNA-guided DNA integration. *Nature*, 571(7764), 219–225.
- Strecker, J. et al. (2019). RNA-guided DNA insertion with CRISPR-associated transposases. *Science*, 365(6448), 48–53.

4. Additional references:

- ° ____
- ° ____
- ° ____

12. REVISION HISTORY

Version	Date	Author	Changes
1.0	____	____	Initial version

13. APPENDICES

Appendix A: Primer Design Worksheet

Project: ____

Target Genomic Locus: ____

Note: gRNA design and cloning is covered in the Gibson Assembly SOP. This worksheet covers only PCR verification primers.

A.1 PCR Screening Primers

External-External Primers (for clonal screening):

Primer Name	Sequence (5' → 3')	Length	Tm	Location
Upstream_Fwd		__ bp	__ °C	~500 bp upstream of target
Downstream_Rev		__ bp	__ °C	~500 bp downstream of integration

Expected Products: - Wild-type (no integration): ~__ bp - Integration: ~__ bp (includes payload)

A.2 Junction Sequencing Primers

Genome-Transposon Junction Primers:

Primer Name	Sequence (5' → 3')	Length	Tm	Target
Genome_Fwd		__ bp	__ °C	Genome, ~300 bp upstream
TnsR_Rev	5'-CGCCAGGGTTTTCCCAGTC-3'	19 bp	60°C	Transposon R end (constant)
TnsL_Fwd	5'-CTGTAGCGCGTTTTTCATCGG-3'	20 bp	60°C	Transposon L end (constant)
Genome_Rev		__ bp	__ °C	Genome, ~300 bp downstream

Expected Products: - 5' junction (GenomeFwd + TnsRRev): ~__ bp - 3' junction (TnsLFwd + GenomeRev): ~__ bp

A.4 Sequencing Primers**For Plasmid Verification:**

Primer Name	Sequence (5' → 3')	Use
CASTarrayseq	5'-GAAACCGTTCGTCATCATCGTC-3'	Verify gRNA spacer in CRISPR array

Appendix B: Reagent Preparation Recipes

All solutions should be filter-sterilized (0.22 µm) or autoclaved as indicated.

Reagent	Recipe	Storage
50% Glycerol (sterile)	Mix 250 mL glycerol + 250 mL MQ water. Autoclave.	4°C, <6 months
10% Glycerol (sterile)	Mix 50 mL glycerol + 450 mL MQ water. Autoclave or filter-sterilize.	4°C, <6 months
1× PBS (pH 7.4)	8 g NaCl, 0.2 g KCl, 1.44 g Na ₂ HPO ₄ , 0.24 g KH ₂ PO ₄ per 1 L water. Autoclave.	RT, 12 months
Spectinomycin (100 mg/mL)	Dissolve 5 g powder in 50 mL water. Filter-sterilize (0.22 µm). Aliquot 1 mL.	-20°C, <12 months
Gentamycin (10 mg/mL)	Dissolve 0.5 g powder in 50 mL water. Filter-sterilize. Aliquot 1 mL.	-20°C, <12 months
DAP (50 mg/mL)	Dissolve 2.5 g diaminopimelic acid in 50 mL water. Filter-sterilize. Aliquot 1 mL.	-20°C, <12 months
TAE Buffer (50×)	242 g Tris base, 57.1 mL glacial acetic acid, 100 mL 0.5 M EDTA (pH 8.0). Adjust to 1 L.	RT, 12 months
1% Agarose Gel	Dissolve 1 g agarose in 100 mL 1× TAE. Microwave until clear. Cool to 60°C, add 10 µL SYBR Safe. Pour.	Use fresh
LB Medium	10 g tryptone, 5 g yeast extract, 10 g NaCl per 1 L water. Autoclave.	RT, 1 month
LB Agar	LB medium + 15 g agar per 1 L. Autoclave. Cool to 55°C, add antibiotics. Pour plates.	4°C, 1 month

Working Antibiotic Concentrations:

Antibiotic	Stock	Working Concentration	Use
Spectinomycin	100 mg/mL	100 µg/mL	Donor selection (for plasmid maintenance)
Selection antibiotic	As appropriate	Determine MIC for recipient	Transconjugant selection
DAP	50 mg/mL	0.3 mM (300 µM)	Donor growth supplement

Note: Determine the minimum inhibitory concentration (MIC) of the selection antibiotic for your recipient strain before starting conjugation experiments.

Appendix C: Data Recording Templates

Template 1: Conjugation Record

Date: __

Donor strain:

Recipient strain:

Donor:Recipient ratio: : __

Mating time: hours at 30°C

Selection: LB + __ (__ µg/mL), no DAP

Colony counts: - Transconjugants (LB + Gent): __ CFU/mL - **Total cells (LB only):** CFU/mL -

Conjugation efficiency: __ %

Result: ☐ Success ☐ Failed

Notes: __

Template 2: PCR Screening Record

Date: __

Primer set: UpstreamFwd + DownstreamRev

Template: Colony lysates from transconjugants

Polymerase: Q5 or OneTaq

Annealing temp: °C

Extension time: __ sec

Cycles: 30

Gel Results: - Expected WT band: __ bp - Expected integration band: bp - Colonies screened: -
Integration-positive colonies: __

Result: ☐ Success ☐ Optimize

Notes: __

Template 3: Sequencing Verification Record

Date: __

Clone ID: Clone #

Junction: ☐ 5' (upstream) ☐ 3' (downstream)

Primer used:

Sequencing service: __

Results: - 5-bp TSD identified: ☐ Yes ☐ No → Sequence: 5'-__-3' - **Distance from target:** __ bp
(expected: ~49 bp) - Orientation: ☐ T-RL ☐ T-LR - Junction sequence correct: ☐ Yes ☐ No

Overall: ☐ PASS ☐ FAIL

Notes: __

Template 4: Glycerol Stock Record

Date: __

Strain name:

Clone ID: Clone #

Parent strain:

Modification: at __

Storage: - Vial 1: Box __, Position __, Freezer - Vial 2: Box __, Position __, Freezer - Vial 3: Box __, Position __, Freezer __

Lab notebook page: _

Electronic files: _

Appendix D: Quick Reference - CAST Integration Mechanism

CAST (CRISPR-Associated Transposase) Integration Features:

1. Target Recognition:

- gRNA: 32-bp target sequence
- PAM: 5'-CN-3' (C = C/G/T/A, N = any) for most Type I-F CAST systems
- TniQ-Cascade complex binds target

2. Integration Site:

- Distance: ~49 bp downstream of target (3' end)
- Variability: ± 5 bp depending on local sequence
- Target site is NOT disrupted

3. Integration Products:

- **T-RL orientation (preferred, >90%):** Right end proximal to target
- **T-LR orientation (rare, <10%):** Left end proximal to target
- **5-bp Target Site Duplication (TSD):** Flanks insertion on both sides

4. Transposon Structure:

- Mini-transposon: Payload flanked by L and R ends
- Payload size: Optimal 0.5-2 kb, functional up to >10 kb
- Transposon ends: ~200 bp (L) and ~60 bp (R, can be truncated to 57 bp)

5. Efficiency:

- Single-target integration: 50-99% (depending on conditions)
- Temperature effect: Lower temp (25-30°C) improves large payloads
- Specificity: >99% on-target for well-designed gRNAs

6. Selection:

- Suicide vector: No plasmid curing needed
 - Integration selection: Payload marker (defined in Project Information)
 - Transconjugants: Selected on antibiotic + no DAP
-

END OF SOP

Performed by: _

Date: _

Signature: _

Reviewed by: _

Date: _

Signature: _